

Gerätezentrum Zytometrie

BD Arial

Filter configuration

Laser	PMT	Filters	Mirrors	Dyes
Blue Laser (488 nm)				
	A	780/60	735LP	PE-Cy7
	B	695/40	655LP	PerCP-Cy5.5 or PI
	C	610/20	595LP	PE-Texas Red
	D	576/26	556LP	PE or PI
	E	530/30	502LP	FITC
	F	488/10	NA	SSC
Red Laser (640 nm)				
	A	780/60	735LP	APC-Cy7
	B	660/20	NA	APC